## CARR® Centritech Separation Systems CARR® Powerfuge® P18



## **Features**

- Drier solids and clear liquids
- Sealed automatic operation
- Full CIP/SIP designs
- cGMP construction
- Designs available for ASME-rated boundary, hazardous locations, and aseptic operations
- · State-of-the-art manufacturing facility located in Clearwater, Florida USA



## CARR®Centritech Separation Systems CARR® Powerfuge® P18

Pneumatic Scale Angelus' unique CARR Powerfuge P18 System provides continuous-flow two-phase separation of solids and liquids at rates up to 1,700 liters/hour for batches up to 10,000 liters.

Using an innovative design, the feed is introduced through a stationary pipe and accelerated to full rotational speed before entering the titanium-alloy bowl. In this system, centrifugal forces as great as 20,000 g permit separation of even sub-micron particles. Compressed solids are periodically removed with a fully-automated scraping cycle. After the solids have been discharged, automated CIP/SIP cycles can be performed via the PLC-controlled operating sequence.

The Powerfuge P18 is designed for ease of operation, with seal feed and drain lines allowing for aerosol-free operation. Similarly, no special tools are required to disassemble the system. Integrated turnkey systems incorporating feed delivery options, CIP and SIP skids, temperature control, and closed-loop centrate pumping under vacuum are available.

The Powerfuge P18 system meets cGMP and CE requirements and has been validated to meet chemical and pharmaceutical industry protocols. The Powerfuge P18 is engineered to optimize ease of operation while maximizing throughput; each step in the process can be readily automated so that time spent on routine maintenance can be minimized.

Optional Service 23 (584) Spindle Cooling High Speed 1" Tri-Clamps: Bearings Housing Wash Feed In (1 1/2") Bowl Case Cooling Out 96 (2439) Cooling In Wash Bottle Wash Centrate Port (4" Tri-Clamp) Solids Port Pneumatic Seal Residual Optional Solids Liquid 40 (1016) 39 (991) 60 (1524)

Note: Dimensions are displayed in inches (millimeters).

For more information, please contact your CARR Centritech representative at:

5320 140th Avenue North Clearwater, FL 33760 USA (727) 535-4100

Otto-Hahn-Strasse 54 D-63303 Dreieich-Sprendlingen Germany (+49) 6103-38090

Angelus | CARR Centritech | Consolidated | Continental/Closetech | Hema Mateer Burt | Multi-Tech Systems | OnTrack | Zepf Change Parts



<sup>\*</sup>Vessels designed for integration with the P18 for vacuum operation and solids handling are available.

<sup>\*\*</sup>Actual capacity (throughput) is dependent on the characteristics of the fluids, amount and type of the solids and required degree of clarification. Practical capacities may differ substantially for your applications.